

## Tip Compatibility Chart

Expell Tip Number	10SF	10LF	20F	100F	200F	1000F
<b>Capp manual pipettes</b>						
	0.2-2ul	0.2-2ul	2-20ul	5-50ul	2-20ul	100-1000ul
	0.5-10ul	0.5-10ul		10-100ul	5-50ul	
					10-100ul	
					20-200ul	
<b>Capp Electronic Pipettes</b>						
			2-20ul		10-200ul	100-1200ul
<b>Brand Transferpette,</b>						
Transferpette S and	0.5-10ul	0.5-10ul	2-20ul		20-200ul	100-1000ul
Transferpette electronic						
<b>Eppendorf Reference,</b>						
Research Plus	0.5-10ul	0.1-2.5ul 0.5-10ul	2-20ul		20-200ul	100-1000ul
<b>Finnpipette Digital</b>						
	0.5-10ul	0.5-10ul	2-20ul	5-50ul	5-50ul 20-200ul	100-1000ul
<b>Gilson Pipetman</b>						
	0.1-2ul	0.1-2ul	2-20ul	2-20ul	2-20ul	100-1000ul
	0.5-10ul	0.5-10ul		10-100ul	10-100ul 20-200ul	
<b>Rainin Classic (non-LTS)</b>						
	0.1-2ul	0.1-2ul	2-20ul	2-20ul	2-20ul	100-1000ul
	0.5-10ul	0.5-10ul		10-100ul	10-100ul 20-200ul	
<b>VWR, Labnet,</b>						
HTL Manual Pipettes	0.1-2ul	0.1-2ul	2-20ul	2-20ul	2-20ul	100-1000ul
	0.5-10ul	0.5-10ul		10-100ul	10-100ul 20-200ul	

BioVentures, Inc,

Website: [www.bioventures.com](http://www.bioventures.com) / email: [support@bioventures.com](mailto:support@bioventures.com)

Phone: 1-800-235-8938/ Fax: 615-896-4837